


Adapting to Olympus Nonvideo, 100 and 130 Series Immersible Endoscopes**



*Part #16030 is a filter insert that must be replaced every 6 to 8 weeks (or more often if required).

CAUTION: Remove all valves (air/water, suction, CO₂) and biopsy caps, selector mechanism, etc.) consult scope manual for proper removal and cleaning procedures of components. Install certain reprocessing parts supplied by Olympus and all adapters supplied at time of In-Service. To confirm proper hook up refer to In-Service Program Report Form and/or Scope Adaptation For All Model Endoscopes Chart. **If you have any questions on adapting to your scopes or if you acquire any new scopes not listed on your In-Service Program/Report Form please call Custom Ultrasonics.**

** This adapter chart for Olympus Flexible Endoscopes illustrates various models and does not represent an actual model and is for reference only. Note: Use *only* Olympus Adapter MH-974 (it has one green loop) or MB-1 with 16008 or 16008M for Specialty Channels.