



# 1 **LOAD** Endoscope

- Upon completion of manual cleaning, place endoscope in basin with control head in right rear corner.

# 2 **ATTACH** Endoscope Hook-up and Leak Tester

- Connect appropriate hook-up and leak tester to basin and endoscope.
- Verify tubing is not kinked or pinched.
- Place floating lid on basin and close reprocessor lid.



# 3 **SELECT** Station

- Press “Station Select” to choose station A or B.

# 4 **ENTER** ID Data

- Press “ID Data” and use the number key pad to enter scope, operator, patient, and physician data (if applicable).
- Press “Enter”.

# 5 **START** Reprocessing Cycle

- Press “Start” to begin reprocessing Cycle.



# 6 **TEST** Disinfectant Concentration

- At the end of the cycle, when prompted, test the MRC with Rapicide® PA test strip.
- Press the HLD Pass or HLD Fail button.

# 7 **REMOVE** Endoscope

- LCD Screen will state “Cycle Complete”.
- Disconnect hook-up connectors and leak tester and remove endoscope.

This Quick Start Guide supplements, but does not replace the information provided in the DSD EDGE™ User Manual, which came with your reprocessor.